



ArchAlignPro

**Your Trusted Partner in Premium
Orthodontic Supplies**



OEM/ODM Services





Advanced Production Equipment



High-standard Production Processes

ArchAlignPro

PRODUCT PORTFOLIO OVERVIEW



Advanced Bracket Systems

INNOVATIVE ORTHODONTIC SOLUTIONS

Traditional Metal Brackets

- Prescription: MBT / Roth
- Type: Standard Metal / Mini Metal
- Material: Medical-grade Stainless Steel
- Models: No Hook, 3-Hooks, 345-Hooks
- Slot Size: 0.018 , 0.022
- Base Type: 80-Mesh Mesh Base



Passive Self-Ligating Brackets

- Prescription: MBT / Roth
- Mechanism: Slide-door Passive Clip
- Material: Medical-grade Stainless Steel
- Slot Size: 0.018 , 0.022
- Feature: Low-friction Sliding Design
- Models: No Hook, 3-Hooks, 345-Hooks



Aesthetic Ceramic Brackets

- Prescription: MBT / Roth
- Material: Ceramic
- Slot Size: 0.018 , 0.022
- Models: No Hook, 3-Hooks, 345-Hooks



Comprehensive Archwire Portfolio

Super-Elastic NiTi Wires

- Features: Excellent shape memory, constant moderate forces, high flexibility.
- Clinical Phase: Initial leveling and aligning phase



Thermal Activated NiTi Wires

- Features: Soft at room temperature, fully activated by intraoral heat, gentle and continuous forces.
- Clinical Phase: Early leveling phase for severe crowding or sensitive patients



Copper NiTi Wires

- Features: Added copper for precise thermal responsiveness, enhanced unloading forces, superior resistance to permanent deformation.
- Clinical Phase: Advanced alignment and early space closure phase



Stainless Steel Wires

- Features: High stiffness, excellent formability for custom bending, low friction.
- Clinical Phase: Mid-to-late stage: Space closure, arch expansion, and torque control



Reverse Curve NiTi Wires

- Features: Pre-formed "rocking chair" curve, exerts vertical intrusive/extrusive forces without complex bending.
- Clinical Phase: Mid-stage: Deep bite correction and leveling of the Spee curve



Orthodontic mini-screw kit

Orthodontic Mini-Implant Screws (TADs)

- **Function:** Provides stable, rigid anchorage in the alveolar bone, enabling precise tooth movement without reliance on neighboring teeth.
- **Ideal Scenarios:** Ideal for complex treatments, including anterior en-masse retraction, molar intrusion, space closure, and occlusal plane correction.



Orthodontic Mini-Implant Driver


- **Function:** Designed for secure, controlled insertion of mini-implants, ensuring optimal torque and tactile feedback for clinical accuracy.
- **Ideal Scenarios:** Essential for all chairside orthodontic TAD placement procedures, offering reliable, efficient handling in various intraoral sites.



Precision Crimpable Hooks



Crimpable Ball Hooks



Crimpable Hooks
(Curved/Double/Triple)



Bondable Lingual Cleat
(with double wings)



Bondable Lingual Button
(with Hook)



Crimpable Stops



Molar Band

Comprehensive Lingual Button Series



Bondable Lingual Buttons
(Rectangular Base)



Round Base Lingual Buttons



Bondable Coil Eyelet



Bondable Lingual Ring



Perforated Lingual Buttons



Eyelet with Chain

Precision Auxiliaries & Components



Bite Turbos/
Tongue Spurs



Five-Loop Torque Springs



Torque Springs



Bondable Buccal Tube



Open Coil Spring/
Closed Coil Spring



Preformed Ligature Ties

Elastics & Clinical Essentials



Elastic Ligature Ties



Elastomeric Chain



Separators



Rotation Wedges



Orthodontic Elastics



Orthodontic Wax

ArchAlignPro

Orthodontic Systems

Shaping Smiles with Precision and Innovation.



Office Address:

Shenzhen, Guangdong Province, China



Global Sales Email: sales@archalignpro.com



Official Website: www.archalignpro.com



Phone/WhatsApp: +853 8684 8787